

Lysosome-associated membrane glycoprotein 1

Catalog No: #AP85032

Package Size: #AP85032-1 50ug #AP85032-2 100ug #AP85032-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosome-associated membrane glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LAMP-1,120 kDa lysosomal membrane glycoprotein,CD107 antigen-like family member A,LGP-120,Lysosomal membrane glycoprotein A,LGP-A,P2B,CD107a,Lamp-1
Accession No.	Uniprot ID:P11438
Uniprot	P11438
GeneID	16783;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamp1

Note: This product is for in vitro research use only