

# Neural cell adhesion molecule L1

Catalog No: #AP85035



Package Size: #AP85035-1 50ug #AP85035-2 100ug #AP85035-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neural cell adhesion molecule L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N-CAM-L1,CD171,CAML1,MIC5
Accession No.	Uniprot ID:P32004
Uniprot	P32004
GeneID	3897;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:L1CAM

Note: This product is for in vitro research use only