

Protein ERGIC-53

Catalog No: #AP85036



Package Size: #AP85036-1 50ug #AP85036-2 100ug #AP85036-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein ERGIC-53
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ER-Golgi intermediate compartment 53 kDa protein, Lectin mannose-binding 1,p58,Ergic53
Accession No.	Uniprot ID:Q62902
Uniprot	Q62902
GeneID	116666;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lman1

Note: This product is for in vitro research use only