

Lumican

Catalog No: #AP85038



Package Size: #AP85038-1 50ug #AP85038-2 100ug #AP85038-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lumican
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Corneal keratan sulfate proteoglycan 37B core protein,Keratan sulfate proteoglycan,KSPG,LDC
Accession No.	Uniprot ID:Q05443
Uniprot	Q05443
GeneID	280847;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LUM

Note: This product is for in vitro research use only