

Leukemia inhibitory factor

Catalog No: #AP85039



Package Size: #AP85039-1 50ug #AP85039-2 100ug #AP85039-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukemia inhibitory factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LIF,Differentiation-stimulating factor,D factor
Accession No.	Uniprot ID:P09056
Uniprot	P09056
GeneID	16878;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lif

Note: This product is for in vitro research use only