

## Laminin subunit alpha-5

Catalog No: #AP85040



Package Size: #AP85040-1 50ug #AP85040-2 100ug #AP85040-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Laminin subunit alpha-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Laminin-10 subunit alpha,Laminin-11 subunit alpha,Laminin-15 subunit alpha,KIAA0533,KIAA1907
Accession No.	Uniprot ID:O15230
Uniprot	O15230
GeneID	3911;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LAMA5

Note: This product is for in vitro research use only