

L-lactate dehydrogenase B chain

Catalog No: #AP85042



Package Size: #AP85042-1 50ug #AP85042-2 100ug #AP85042-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	LDH-B,LDH heart subunit,LDH-H
Accession No.	Uniprot ID:P00336
Uniprot	P00336
GeneID	100621540;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LDHB

Note: This product is for in vitro research use only