

Hepatic triacylglycerol lipase

Catalog No: #AP85043



Package Size: #AP85043-1 50ug #AP85043-2 100ug #AP85043-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatic triacylglycerol lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HL,Lipase member C
Accession No.	Uniprot ID:Q3SZ79
Uniprot	Q3SZ79
GeneID	535192;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LIPC

Note: This product is for in vitro research use only