

Creatine kinase B-type

Catalog No: #AP85044



Package Size: #AP85044-1 50ug #AP85044-2 100ug #AP85044-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Creatine kinase B-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B-CK,Creatine kinase B chain,CKBB
Accession No.	Uniprot ID:P12277
Uniprot	P12277
GeneID	1152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CKB

Note: This product is for in vitro research use only