

Ribosomal protein S6 kinase beta-1

Catalog No: #AP85045



Package Size: #AP85045-1 50ug #AP85045-2 100ug #AP85045-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal protein S6 kinase beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	S6K-beta-1,70 kDa ribosomal protein S6 kinase 1,P70S6K1,Ribosomal protein S6 kinase I,Serine/threonine-protein kinase 14A,p70 ribosomal S6 kinase alpha,p70 S6 kinase alpha,STK14A
Accession No.	Uniprot ID:P23443
Uniprot	P23443
GeneID	6198;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RPS6KB1

Note: This product is for in vitro research use only