

Lupus La protein homolog

Catalog No: #AP85046



Package Size: #AP85046-1 50ug #AP85046-2 100ug #AP85046-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Lupus La protein homolog |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | La autoantigen homolog,La ribonucleoprotein,Ss-b |
| Accession No. | Uniprot ID:P38656 |
| Uniprot | P38656 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ssb

Note: This product is for in vitro research use only