

# La-related protein 7

Catalog No: #AP85051



Package Size: #AP85051-1 50ug #AP85051-2 100ug #AP85051-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	La-related protein 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	La ribonucleoprotein domain family member 7,P-TEFb-interaction protein for 7SK stability,PIP7S,HDCMA18P
Accession No.	Uniprot ID:Q4G0J3
Uniprot	Q4G0J3
GeneID	51574;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LARP7

Note: This product is for in vitro research use only