

Galectin-3-binding protein

Catalog No: #AP85053



Package Size: #AP85053-1 50ug #AP85053-2 100ug #AP85053-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Lectin galactoside-binding soluble 3-binding protein
Accession No.	Uniprot ID:A7E3W2
Uniprot	A7E3W2
GeneID	531137;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LGALS3BP

Note: This product is for in vitro research use only