

Interleukin-21

Catalog No: #AP85054



Package Size: #AP85054-1 50ug #AP85054-2 100ug #AP85054-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IL-21
Accession No.	Uniprot ID:Q76LU5
Uniprot	Q76LU5
GeneID	378475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL21

Note: This product is for in vitro research use only