

Melanoma-associated antigen 4

Catalog No: #AP85061



Package Size: #AP85061-1 50ug #AP85061-2 100ug #AP85061-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanoma-associated antigen 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cancer/testis antigen 1.4,CT1.4,MAGE-4 antigen,MAGE-41 antigen,MAGE-X2 antigen,MAGE4
Accession No.	Uniprot ID:P43358
Uniprot	P43358
GeneID	4103;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAGEA4

Note: This product is for in vitro research use only