

Induced myeloid leukemia cell differentiation protein Mcl-1 homolog



Catalog No: #AP85063

Orders: order@signalwayantibody.com

Package Size: #AP85063-1 50ug #AP85063-2 100ug #AP85063-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Induced myeloid leukemia cell differentiation protein Mcl-1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q9Z1P3
Uniprot	Q9Z1P3
GeneID	60430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mcl1

Note: This product is for in vitro research use only