

UMP-CMP kinase

Catalog No: #AP85066



Package Size: #AP85066-1 50ug #AP85066-2 100ug #AP85066-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | UMP-CMP kinase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Deoxycytidylate kinase,CK,Nucleoside-diphosphate kinase,Uridine monophosphate/cytidine monophosphate kinase,UMP/CMP kinase,CMPK |
| Accession No. | Uniprot ID:Q2KIW9 |
| Uniprot | Q2KIW9 |
| GeneID | 509965; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CMK1

Note: This product is for in vitro research use only