

Interleukin-12 subunit alpha

Catalog No: #AP85072



Package Size: #AP85072-1 50ug #AP85072-2 100ug #AP85072-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-12 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-12A,Cytotoxic lymphocyte maturation factor 35 kDa subunit,CLMF p35,IL-12 subunit p35,NK cell stimulatory factor chain 1,NKSF1,NKSF1
Accession No.	Uniprot ID:P29459
Uniprot	P29459
GeneID	3592;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL12A

Note: This product is for in vitro research use only