

Interleukin-6 receptor subunit beta

Catalog No: #AP85075



Package Size: #AP85075-1 50ug #AP85075-2 100ug #AP85075-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-6 receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-6 receptor subunit beta, Interleukin-6 signal transducer, Membrane glycoprotein 130, gp130, Oncostatin-M receptor subunit alpha, CD130
Accession No.	Uniprot ID:Q00560
Uniprot	Q00560
GeneID	16195;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL6st

Note: This product is for in vitro research use only