

## Kininogen-1

Catalog No: #AP85076



Package Size: #AP85076-1 50ug #AP85076-2 100ug #AP85076-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kininogen-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-2-thiol proteinase inhibitor,Fitzgerald factor,High molecular weight kininogen,HMWK,Williams-Fitzgerald-Flaujeac factor,BDK,KNG
Accession No.	Uniprot ID:P01042
Uniprot	P01042
GeneID	3827;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:KNG1

Note: This product is for in vitro research use only