

Low-density lipoprotein receptor

Catalog No: #AP85082



Package Size: #AP85082-1 50ug #AP85082-2 100ug #AP85082-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	LDL receptor
Accession No.	Uniprot ID:Q28832
Uniprot	Q28832
GeneID	396801;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LDLR

Note: This product is for in vitro research use only