

Myelin and lymphocyte protein

Catalog No: #AP85084



Package Size: #AP85084-1 50ug #AP85084-2 100ug #AP85084-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myelin and lymphocyte protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	T-lymphocyte maturation-associated protein
Accession No.	Uniprot ID:O09198
Uniprot	O09198
GeneID	17153;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mal

Note: This product is for in vitro research use only