

NADP-dependent malic enzyme, mitochondrial

Catalog No: #AP85085



Package Size: #AP85085-1 50ug #AP85085-2 100ug #AP85085-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADP-dependent malic enzyme, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NADP-ME, Malic enzyme 3
Accession No.	Uniprot ID:Q8BMF3
Uniprot	Q8BMF3
GeneID	109264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Me3

Note: This product is for in vitro research use only