

Latent-transforming growth factor beta-binding protein

2

Catalog No: #AP85087

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85087-1 50ug #AP85087-2 100ug #AP85087-3 1mg

Description

Product Name	Latent-transforming growth factor beta-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LTBP-2,C14orf141,LTBP3
Accession No.	Uniprot ID:Q14767
Uniprot	Q14767
GeneID	4053;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LTBP2

Note: This product is for in vitro research use only