

Ly6/PLAUR domain-containing protein 3

Catalog No: #AP85088



Package Size: #AP85088-1 50ug #AP85088-2 100ug #AP85088-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Ly6/PLAUR domain-containing protein 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | GPI-anchored metastasis-associated protein C4.4A,C4.4a |
| Accession No. | Uniprot ID:O55162 |
| Uniprot | O55162 |
| GeneID | 60378; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lypd3

Note: This product is for in vitro research use only