

## Keratin, type I cytoskeletal 18

Catalog No: #AP85089



Package Size: #AP85089-1 50ug #AP85089-2 100ug #AP85089-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type I cytoskeletal 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytokeratin endo B, Keratin D, Cytokeratin-18, CK-18, Keratin-18, K18, Kerd, Krt1-18
Accession No.	Uniprot ID:P05784
Uniprot	P05784
GeneID	16668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Krt18

Note: This product is for in vitro research use only