

Mannan-binding lectin serine protease 2

Catalog No: #AP85092



Package Size: #AP85092-1 50ug #AP85092-2 100ug #AP85092-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mannan-binding lectin serine protease 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MBL-associated serine protease 2, Mannose-binding protein-associated serine protease 2, MASP-2
Accession No.	Uniprot ID:Q9JJS8
Uniprot	Q9JJS8
GeneID	64459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Masp2

Note: This product is for in vitro research use only