

Mitotic spindle assembly checkpoint protein MAD1

Catalog No: #AP85093



Package Size: #AP85093-1 50ug #AP85093-2 100ug #AP85093-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitotic spindle assembly checkpoint protein MAD1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitotic arrest deficient 1-like protein 1,MAD1-like protein 1, Mad1
Accession No.	Uniprot ID:Q9WTX8
Uniprot	Q9WTX8
GeneID	17120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mad111

Note: This product is for in vitro research use only