

## Lysophosphatidic acid receptor 4

Catalog No: #AP85095



Package Size: #AP85095-1 50ug #AP85095-2 100ug #AP85095-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysophosphatidic acid receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LPA receptor 4,G-protein coupled receptor 23,P2Y purinoceptor 9,P2Y9,Purinergic receptor 9,Gpr23,Lpa4,P2y9
Accession No.	Uniprot ID:Q8BLG2
Uniprot	Q8BLG2
GeneID	78134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lpar4

Note: This product is for in vitro research use only