

# Lutropin subunit beta

Catalog No: #AP85096



Package Size: #AP85096-1 50ug #AP85096-2 100ug #AP85096-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lutropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Luteinizing hormone subunit beta,LH-B,Lutropin beta chain
Accession No.	Uniprot ID:P04651
Uniprot	P04651
GeneID	280839;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LHB

Note: This product is for in vitro research use only