

Leukotriene C4 synthase

Catalog No: #AP85097



Package Size: #AP85097-1 50ug #AP85097-2 100ug #AP85097-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukotriene C4 synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LTC4 synthase,Leukotriene-C(4) synthase
Accession No.	Uniprot ID:Q2NKS0
Uniprot	Q2NKS0
GeneID	511749;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LTC4S

Note: This product is for in vitro research use only