

Lymphocyte antigen 6C1

Catalog No: #AP85098



Package Size: #AP85098-1 50ug #AP85098-2 100ug #AP85098-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphocyte antigen 6C1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ly-6C1
Accession No.	Uniprot ID:P0CW02
Uniprot	P0CW02
GeneID	17067;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ly6c1

Note: This product is for in vitro research use only