

## L-selectin

Catalog No: #AP85100



Package Size: #AP85100-1 50ug #AP85100-2 100ug #AP85100-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD62 antigen-like family member L,Leukocyte adhesion molecule 1,LAM-1,Leukocyte-endothelial cell adhesion molecule 1,LECAM1,Lymph node homing receptor,Lymphocyte antigen 22,Ly-22,Lymphocyte surface MEL-14 antigen,CD62L,Lnhr,Ly-22
Accession No.	Uniprot ID:P30836
Uniprot	P30836
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sell

Note: This product is for in vitro research use only