

Laminin subunit alpha-1

Catalog No: #AP85104



Package Size: #AP85104-1 50ug #AP85104-2 100ug #AP85104-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Laminin subunit alpha-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Laminin A chain,Laminin-1 subunit alpha,Laminin-3 subunit alpha,S-laminin subunit alpha,S-LAM alpha,Lama,Lama-1 |
| Accession No. | Uniprot ID:P19137 |
| Uniprot | P19137 |
| GeneID | 16772; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lama1

Note: This product is for in vitro research use only