

Protein LYRIC

Catalog No: #AP85110



Package Size: #AP85110-1 50ug #AP85110-2 100ug #AP85110-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein LYRIC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3D3/LYRIC,Astrocyte elevated gene-1 protein,AEG-1,Lysine-rich CEACAM1 co-isolated protein,Metadherin,Metastasis adhesion protein,AEG1,LYRIC
Accession No.	Uniprot ID:Q86UE4
Uniprot	Q86UE4
GeneID	92140;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MTDH

Note: This product is for in vitro research use only