

Mitogen-activated protein kinase kinase kinase 7

Catalog No: #AP85111



Package Size: #AP85111-1 50ug #AP85111-2 100ug #AP85111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase kinase kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Transforming growth factor-beta-activated kinase 1,TGF-beta-activated kinase 1,Tak1
Accession No.	Uniprot ID:Q62073
Uniprot	Q62073
GeneID	26409;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map3k7

Note: This product is for in vitro research use only