

Multidrug resistance protein 1

Catalog No: #AP85115



Package Size: #AP85115-1 50ug #AP85115-2 100ug #AP85115-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATP-binding cassette sub-family B member 1,P-glycoprotein 1,CD243,Mdr1,Mdr1b,Pgy1
Accession No.	Uniprot ID:P43245
Uniprot	P43245
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcb1

Note: This product is for in vitro research use only