

Multidrug resistance protein 2

Catalog No: #AP85116



Package Size: #AP85116-1 50ug #AP85116-2 100ug #AP85116-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-binding cassette sub-family B member 4,P-glycoprotein 2,Mdr2,Pgy-2,Pgy2
Accession No.	Uniprot ID:P21440
Uniprot	P21440
GeneID	18670;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcb4

Note: This product is for in vitro research use only