

Late cornified envelope protein 3D

Catalog No: #AP85119



Package Size: #AP85119-1 50ug #AP85119-2 100ug #AP85119-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Late cornified envelope protein 3D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Late envelope protein 16, Small proline-rich-like epidermal differentiation complex protein 6A, Small proline-rich-like epidermal differentiation complex protein 6B, LEP16, SPRL6A, SPRL6B
Accession No.	Uniprot ID:Q9BYE3
Uniprot	Q9BYE3
GeneID	84648;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LCE3D

Note: This product is for in vitro research use only