

Krueppel-like factor 5

Catalog No: #AP85120



Package Size: #AP85120-1 50ug #AP85120-2 100ug #AP85120-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Krueppel-like factor 5 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Basic transcription element-binding protein 2,BTE-binding protein 2,Colon krueppel-like factor,GC-box-binding protein 2,Intestinal-enriched krueppel-like factor,Transcription factor BTEB2,BTEB2,CKLF,IKLF |
| Accession No. | Uniprot ID:Q13887 |
| Uniprot | Q13887 |
| GeneID | 688; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KLF5

Note: This product is for in vitro research use only