

Lumican

Catalog No: #AP85125



Package Size: #AP85125-1 50ug #AP85125-2 100ug #AP85125-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lumican
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Keratan sulfate proteoglycan lumican,KSPG lumican,LDC,SLRR2D
Accession No.	Uniprot ID:P51884
Uniprot	P51884
GeneID	4060;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LUM

Note: This product is for in vitro research use only