

## P-selectin

Catalog No: #AP85126



Package Size: #AP85126-1 50ug #AP85126-2 100ug #AP85126-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	P-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD62 antigen-like family member P,Granule membrane protein 140,GMP-140,Leukocyte-endothelial cell adhesion molecule 3,LECAM3,Platelet activation dependent granule-external membrane protein,PADGEM,CD62P,Gmp
Accession No.	Uniprot ID:Q01102
Uniprot	Q01102
GeneID	20344;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Selp

Note: This product is for in vitro research use only