

Gastric triacylglycerol lipase

Catalog No: #AP85127



Package Size: #AP85127-1 50ug #AP85127-2 100ug #AP85127-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gastric triacylglycerol lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GL
Accession No.	Uniprot ID:P07098
Uniprot	P07098
GeneID	8513;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LIPF

Note: This product is for in vitro research use only