

# Lysyl oxidase homolog 1

Catalog No: #AP85129



Package Size: #AP85129-1 50ug #AP85129-2 100ug #AP85129-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysyl oxidase homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lysyl oxidase 2,Lysyl oxidase-like protein 1,Lox2,Loxl
Accession No.	Uniprot ID:P97873
Uniprot	P97873
GeneID	16949;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Loxl1

Note: This product is for in vitro research use only