

Leucine-rich repeat and calponin homology domain-containing protein 4

Catalog No: #AP85130

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85130-1 50ug #AP85130-2 100ug #AP85130-3 1mg

Description

Product Name	Leucine-rich repeat and calponin homology domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q921G6
Uniprot	Q921G6
GeneID	231798;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lrch4

Note: This product is for in vitro research use only