

Transcription factor 4

Catalog No: #AP85131



Package Size: #AP85131-1 50ug #AP85131-2 100ug #AP85131-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	TCF-4,Immunoglobulin transcription factor 2,ITF-2,Thyroglobulin promoter transcription factor TFE,ITF2,TFE
Accession No.	Uniprot ID:P15881
Uniprot	P15881
GeneID	403949;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TCF4

Note: This product is for in vitro research use only