

Transcription factor 4

Catalog No: #AP85132



Package Size: #AP85132-1 50ug #AP85132-2 100ug #AP85132-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	TCF-4, Class A helix-loop-helix transcription factor GE2
Accession No.	Uniprot ID:Q90683
Uniprot	Q90683
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TCF4

Note: This product is for in vitro research use only