

Inter-alpha-trypsin inhibitor heavy chain H3

Catalog No: #AP85133



Package Size: #AP85133-1 50ug #AP85133-2 100ug #AP85133-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inter-alpha-trypsin inhibitor heavy chain H3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ITI heavy chain H3
Accession No.	Uniprot ID:Q61704
Uniprot	Q61704
GeneID	16426;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Itih3

Note: This product is for in vitro research use only