

Malignant T-cell-amplified sequence 1

Catalog No: #AP85136



Package Size: #AP85136-1 50ug #AP85136-2 100ug #AP85136-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Malignant T-cell-amplified sequence 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MCT-1, Multiple copies T-cell malignancies 1
Accession No.	Uniprot ID:Q9DB27
Uniprot	Q9DB27
GeneID	68995;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mcts1

Note: This product is for in vitro research use only