

Galectin-3

Catalog No: #AP85137



Package Size: #AP85137-1 50ug #AP85137-2 100ug #AP85137-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Galectin-3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Gal-3,35 kDa lectin,Carbohydrate-binding protein 35,CBP 35,Galactose-specific lectin 3,IgE-binding protein,Laminin-binding protein,Lectin L-29,Mac-2 antigen |
| Accession No. | Uniprot ID:P08699 |
| Uniprot | P08699 |
| GeneID | 83781; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lgals3

Note: This product is for in vitro research use only